

TriLok™

MEDIAL UNLOADING AND LATERAL STABILITY

The key to the TriLok's® outstanding effectiveness is the combination of the patented FootLok® Strap and the Figure-8 Stirrup Strap. The straps work together to provide unrivaled lateral stability while allowing a full range of motion through plantar and dorsi-flexion. In addition, the FootLok® Strap may be positioned to unload the medial column. The TriLok® is also amazingly thin and low-profile, fitting easily into almost any shoe.

INDICATIONS

- Ankle Sprains
Prevention and Treatment
- Ankle Instability
- Mechanically Support Flexible Flat Foot
- Plantar Fasciitis
- Posterior Tibial Tendon Dysfunction (PTTD)
- Peroneal Tendon Strains
- Subtalar and Ankle Arthritis



UNEQUALED VERSATILITY

LATERAL STABILIZATION



The TriLok™ is unique in its ability to deal with issues on both the lateral and medial sides of the ankle. When positioned laterally, the FootLok™ Strap supports and protects the lateral ligaments and the peroneal tendons. When wrapped medially, the FootLok™ Strap may be positioned to support the plantar fascia or the posterior tibial tendon. No other brace offers such versatility.

The TriLok is effective with almost any type of footwear or with no footwear at all.

MEDIAL COLUMN UNLOADING



Order Numbers

Shoe Size:	Men's	Women's	Europe	UK	Model No.
Small	5-6.5	6-7.5	38-39	4-5.5	53601
Medium	7-9.5	8-11	40-42	6-8.5	53602
Large	10-12.5		43-45	9-11.5	53603
X-Large	13-16		46-49	12-16	53604





“The TriLok® is by far the best ankle brace I have ever used. I can cut hard and play aggressive and my ankle isn’t going anywhere.”

-Stephen B.

ANKLE BRACING



The AFTR DC™, AFTR™ and AFTR with Gel™ provide support and stabilization for the non-operative and post-operative treatment of many foot and ankle pathologies. Both braces include medial and lateral plastic stays as well as figure-8 stirrup straps to limit movement in the ankle and hindfoot.

Both are made from the Bio Skin® material so they perfectly contour to the foot during normal volume changes following injuries and surgery. Bio Skin® is hypo-allergenic, breathable, and thin.

AFTR DC™ DORSAL CLOSURE

The AFTR DC™ features a dorsal closure design making it suitable for:

- When initial foot volume is high, but expected to decrease significantly
- When foot volume is expected to fluctuate during the day
- Early post-op use when wounds are vulnerable

The AFTR DC™ includes adjustable medial and lateral silicone gel pads that wrap behind the malleoli.

Style	Model No.
AFTR DC™ XSmall/Small	53861
AFTR DC™ Medium/Large	53863
AFTR DC™ XLarge/XXLarge	53865

INDICATIONS

- Peroneal/Post-Tib Tendon Repair
- Lateral Ankle Repair
- Fusions
 - Midfoot
 - Hindfoot
 - Ankle
- Total Ankle Replacement
- Ankle Sprains
- Bunionectomy
- Forefoot Procedures
- Severe Edema

INDICATIONS

- Ankle Sprains
- Lateral Ankle Repair
- Hindfoot Arthritis
- Peroneal/Post-Tib Tendon Repair
- Bunionectomy
- Forefoot Procedures



AFTR™ AFTR with Gel™

The AFTR™ and AFTR with Gel™ features a simple, pull-on design making it suitable for:

- Patients with normal or close-to-normal foot volume
- Patients who prefer a simpler application process


The AFTR with Gel™ includes adjustable medial and lateral silicone gel pads that wrap behind the malleoli.

Style	Model No.
AFTR™ XSmall/Small	53801
AFTR™ Medium/Large	53803
AFTR™ XLarge/XXLarge	53805
AFTR with Gel™ XSmall/Small	53831
AFTR with Gel™ Medium/Large	53833
AFTR with Gel™ XLarge/XXLarge	53835

AFTR™ and AFTR DC™ Sizing

Ankle Measurement (*circumference from heel to top of foot - back to heel*)

	inches	centimeters
XSmall/Small	11-13	27.9-33
Medium/Large	12.5-14	31.8-35.6
XLarge/XXLarge	13.5-15	34.3-38.1



WEIL DIGIT STRAPPING SYSTEM™



HCS™ HALLUX CONTROL STRAP™

Splinting The HCS™ may be used post-bunionectomy to maintain the hallux in a neutral position during critical soft-tissue healing.

Range of Motion After primary healing, the HCS™ may be used day or night to stretch connective tissue and increase range of motion.

Exercising The HCS™ can be positioned to provide resistance for FHL strengthening.

Shown with Midfoot Compression Wrap, also compatible with AFTR DC,™ AFTR™ and AFTR with Gel.™

Package	Model No.
Size Small Single	95391
Size Small 10-Pack	95392
Size Large Single	95393
Size Large 10-Pack	95394



WEIL OSTEOTOMY STRAP™

The Weil Osteotomy Strap™ helps prevent dorsal drift and floating toes by splinting the toe in a neutral or plantar-flexed position. The Weil Osteotomy Strap™ can easily be worn in a walking boot, post-op shoe, or normal footwear and may also be used non-operatively to treat metatarsalgia or plantar plate tears.

Shown with Midfoot Compression Wrap, also compatible with AFTR DC,™ AFTR™ and AFTR with Gel.™

Package	Model No.
Single	95371
10-Pack	95372



WEIL EXERCISER™

The Weil Exerciser™ is applied to the lesser digits to provide a home therapy program to strengthen the plantarflexors.

Each set comes with one high-tension and one low-tension exerciser to help patients regain normal FDL strength.

Shown with Midfoot Compression Wrap, also compatible with AFTR DC,™ AFTR™ and AFTR with Gel.™

Package	Model No.
Low/High Tension Single	95381
Low/High Tension 10-Pack	95382
Low Tension Single	95383
Low Tension 10-Pack	95384
High Tension Single	95385
High Tension 10-Pack	95386



WEIL DIGIT STRAPPING ATTACHMENT OPTIONS



CALIBRATED MIDFOOT COMPRESSION WRAP™

The Calibrated Midfoot Compression Wrap™ allows you to use any of the Weil Digit Straps™ without covering the entire foot. The material is treated with a hypoallergenic non-slip finish on the inside of the material to effectively hold it in place. The low profile is well-suited for use inside normal footwear.

Package	Model No.
Single	95361
10-Pack	95362



AFTR DC™

The AFTR DC™ provides mechanical stability and adjustable compression to support healing after a variety of injuries and surgical procedures. The AFTR DC™ is compatible with any of the Weil Digit Straps™.

Style	Model No.
AFTR DC™ XSmall/Small	53861
AFTR DC™ Medium/Large	53863
AFTR DC™ XLarge/XXLarge	53865



AFTR™ AFTR with Gel™

The AFTR™ and AFTR with Gel™ provide mechanical stability with rigid stirrup straps as well as compression for edema prevention and improved proprioception. The AFTR™ and AFTR with Gel™ are compatible with the Hallux Control Strap™ and the Weil Osteotomy Strap™.

Style	Model No.
AFTR™ XSmall/Small	53801
AFTR™ Medium/Large	53803
AFTR™ XLarge/XXLarge	53805
AFTR with Gel™ XSmall/Small	53831
AFTR with Gel™ Medium/Large	53833
AFTR with Gel™ XLarge/XXLarge	53835