

THE WORLD'S FIRST HYBRID KNEE JUST GOT BETTER:
THE ALL NEW 2018 CROSSOVER HYBRID PROSTHETIC KNEE

2018:
New Integrated Design.

Patented:
highly adjustable Tendon and Two-Position Quick Release X/CAM System.

Walk Mode:
New rubberized Protective Cover.

X/CAMs:
allow the highest tuning of tendon response within further progression of knee flexion.

CrossOver Tendons:
provide a return mechanism that is variable with infinitely tunable resistance.

No Tools:
make quick changes "on the fly"

WALK MODE

SPORT MODE

Crossflex™ wedges: act like knee cartilage. Three-pack pre-flexion foam wedges for sport activities pick the wedge degree to suit the activity.

VED bumpers: 75° and 85° durometer ratings. Interchange bumpers to address knee performance during stance phase.

Durable Bike Tech:
Crossover utilizes a highly adjustable pneumatic-hydraulic shock with user adjusted air pressure.

Infinite Possibilities.

Designed To Accommodate Changes.

CrossOver can be used for everyday walking and -with simple additions/adjustments- can be used for low and high activity sports.

Take A Walk On The Wild Side.

Walking knees are good for walking only and sport knees are not good for walking. That's a big problem for the user with real life transitions from one activity to the next. CrossOver bridges that gap, the world's first hybrid knee that delivers high performance in all activities.



Brian Bartlett inventor and athlete.

Adapts To You. No Tools Required.



TOOL FREE

You will be amazed at what you can do with this knee! CrossOver allows you to make adjustments on the spot with no tools. Shred the hill and then walk to the lodge! CrossOver's got all daily aspects of your life covered.

Demo A CrossOver Knee Today!

You or patient can evaluate this extensible knee for 10 days before committing to a purchase. Apply today. Supplies are limited.



VISIT FABTECHSYSTEMS.COM FOR MORE INFORMATION. OR CALL: 800.FABTECH